

TOWER

MEDICAL SYSTEMS LTD.

917-11 Lincoln Avenue

FIXED C-ARM TABLE MODEL FCA-1000

INEXPENSIVE, STRONG AND STURDY TABLES FOR
FLUOROSCOPIC IMAGING AND PAIN CARE IMAGING PROCEDURES



Our Fixed C-Arm Table is an economical yet strong and durable table designed to be used for fluoroscopic imaging procedures as well as pain care imaging studies. This table is **32 ½" high (82 cm)** from floor to fiber resin table top. (excluding mattress) This table is also available extra high at 36.5". The **24" width (50 cm)** and **extra long 78" length (198 cm)** accommodates greater imaging area. This table is rated at 650 lb. (294 Kg) weight capacity. The metal frame is powder coated in a texture gray color and comes standard with a 1" thick black mattress. An optional IV Pole with clamps can easily be attached.

MODEL FCA-1000 Fixed C-Arm Table.....	\$2,850.00
Model FCA-1000-EX (36.5").....	\$3150.00
Optional IV Pole with Clamps	\$150.00

(631) 699-3200

Tower Medical Systems Ltd.

Fax (631) 699-3201

Visit us on the Web at www.towermedicalsystems.com

Or E-Mail to TowerMedicalLtd@AOL.Com