

TOWER

MEDICAL SYSTEMS LTD.

917-11 Lincoln Avenue
Holbrook, NY 11741

SCANNER CART MODEL SC-500

INEXPENSIVE, STRONG AND STURDY TABLES



Our Scanner Cart is an economical yet strong and durable table designed to be used for imaging procedures as well as pain care imaging studies. This table is **32 ½" high (77.5 cm)** from floor to ¼ inch thick phenolic table top. It has a **24" width (50 cm)** and **78" length (188 cm)** and is rated at **450 lb. (203 Kg) weight capacity**. The metal frame is powder coated in a texture gray color and comes standard with 5" locking casters. Options include a 1" or 2" thick Lectrolite Pad (helps control electrostatic discharge), Patient Guard Rails and an IV Pole which are all easily attached.

| | |
|--------------------------------|------------|
| MODEL SC-500 Scanner Cart..... | \$2,150.00 |
| 1" Thick Lectrolite Pad..... | \$269.00 |
| 2" Thick Lectrolite Pad..... | \$305.00 |
| Side Rails..... | \$300.00 |
| Ankle Stirrups | \$380.00 |

(631) 699-3200

Tower Medical Systems Ltd.

Fax (631) 699-3201

Visit us on the Web at www.towermedicalsystems.com

Or E-Mail to Info@TowerMedicalSystems.Com