

Multi-Purpose Power Exam Table

MD250



917-11 Lincoln Avenue
Holbrook, New York 11741
tel. 631.699.3200
fax. 631.699.3201

e-mail info@towermedicalsystems.com
www.towermedicalsystems.com



MD-250

At Tower Medical Systems, the effectiveness, efficiency, and safety of our tables are our main concern. The quality and structure of our product correlates with both the physicians' and patients' needs. By keeping these needs in mind, the table stays true to its purpose, which is to provide comfort to the patient and accessibility for the physician.

The MD-250 is a multipurpose table designed for Ultrasound procedures including proctology, urology, thyroid and OB/GYN (when stirrups are attached). The sleek size of the MD-250 in the upright position makes it ideal for tight quarters. The MD-250 is the first of its kind all-around mobile procedure table. The U-shape front frame assembly, along with the hideaway leg extension, makes this table perfect for all OB/GYN, prostate, and scanning procedures. The addition of the Drain Bag Hoop Assembly facilitates in urological procedures.

The motor-driven MD-250 is manipulated by handheld remote control or foot controls. Three (3) precision Swiss made motors control height, Trendelenberg, and Fowler back positions. The unique telescopic sheet metal design encloses all motorized parts for the safety of both physician and patient. The table uses standard 110V AC 50/60 Hz. Alterations can be made to meet International Specifications.

Support construction is of Fibersin and Velcro® for easy cleaning and the assurance of no warping. The Fibersin used in our tables is certified to meet all requirements of radiation performance standards of 21 CFR Subchapter J.

Patients are secured with 2-inch hook and loop fastener restraining straps and optional swing down guard rails.

MD-250E Echocardiography Table

Our MD-250E Echocardiography Table has a 7 x 10-1/2 inch patient cutout for unobstructed access to the patient's left thorax area. This table enables sonographers to quickly achieve an optimal imaging window, there by shortening exam time.



**Echo Option
MD-250E**

Available Options

*Detachable
Armrests*



Guardrails



*Ankle
Stirrups
Attached*



*Knee
Crutches*



Wedge



*Ankle
Stirrups*



IV Pole also an available option.

*Upholstery: Available in standard sea foam, light blue, mauve, teal, black, forest green, dark grey, and medium grey.

Specifications:

- Patient Capacity: 450 lbs.
- Retractable Foot Rest
- Standard 110V AC 50/60 Hz with battery back up system
- Table Size: 22"x 74"
- 2" Cushions
- Height Adjustability: 28-37 inches
- 5" Locking Casters
- Trendelenberg Positioning +/- 15 degrees
- Motorized Fowler Back 0-85 degrees
- Hand and Foot Controls Standard



* Both height adjustability and Trendelenberg positioning are controlled via hand and foot control.

Low Lease Rates Available

